

MIER2 antibody



Catalog Number: 112617

Product name

MIER2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MIER2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 48-50kd

Preparation

This antibody was obtained by immunization of MIER2 recombinant protein (Accession Number: NM_017550). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112617(MIER2 antibody) at dilution of 1:500