

# MINK1 antibody

Catalog Number: 112623

## Product name

MINK1 antibody

WB: 1:500-1:5000

IP: 1:200-1:1000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

MINK1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human lymphoma tissue, human skeletal muscle tissue, mouse brain tissue. Observed molecular weight by Western-blot: 150 kDa

## Validations




## Preparation

This antibody was obtained by immunization of MINK1 recombinant protein (Accession Number: NM\_153827). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112623(MINK1 antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:

IP Result of anti-MINK1 (IP:Catalog No:112623, 4ug; Detection:Catalog No:112623 1:300) with mouse brain tissue lysate 4000ug.




# MINK1 antibody

Catalog Number: 112623

---

Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:112623(MINK1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:112623(MINK1 Antibody) at dilution of 1:50 (under 40x lens)