

# MIOS antibody

Catalog Number: 112624

## Product name

MIOS antibody

## Specificity

Human; other species not tested.

## Antibody description

MIOS Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 90-100 kDa

## Preparation

This antibody was obtained by immunization of MIOS recombinant protein (Accession Number: BC005883). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112624(MIOS Antibody) at dilution of 1:600



IP Result of anti-MIOS (IP:Catalog No:112624, 4ug; Detection:Catalog No:112624 1:300) with HeLa cells lysate 4000ug.