

MEKK2 antibody

Catalog Number: 112650

Product name

MEKK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEKK2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, K-562 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 65-75 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006609). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112650(MAP3K2 antibody) at dilution of



1:1000

IP Result of anti-MAP3K2 (IP:Catalog No:112650, 4ug; Detection:Catalog No:112650 1:300) with HEK-293 cells lysate 3200ug.