

MLANA antibody



Catalog Number: 112672

Product name

MLANA antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

MLANA Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human malignant melanoma tissue, human skin tissue.

Preparation

This antibody was obtained by immunization of MLANA recombinant protein (Accession Number: NM_005511). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112672(MLANA Antibody) at dilution of 1:700



Immunohistochemistry of paraffin-embedded human malignant melanoma tissue slide using Catalog No:112672(MLANA Antibody) at dilution of 1:50 (under 40x lens)