

MLN antibody

Catalog Number: 112683

Product name

MLN antibody

Specificity

human; other species not tested.

Antibody description

MLN Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue, COLO 320 cells, human colon tissue, human ileum tissue. Observed molecular weight by Western-blot: 13kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002418). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:112683(MLN antibody) at dilution of 1:4000