

MOBKL2A antibody



Catalog Number: 112721

Product name

MOBKL2A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MOBKL2A Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, Jurkat cells, mouse lymph tissue, Raji cells. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of MOBKL2A recombinant protein (Accession Number: XM_047438169). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:112721(MOBKL2A Antibody) at dilution of 1:1000