

MOGAT2 antibody



Catalog Number: 112730

Product name

MOGAT2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MOGAT2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse colon tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 36-38 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_025098). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112730(MOGAT2 antibody) at dilution of 1:300



IP Result of anti-MOGAT2 (IP:Catalog No:112730, 4ug; Detection:Catalog No:112730 1:300) with HEK-293 cells lysate 4800ug.