

MOSPD3 antibody



Catalog Number: 112737

Product name

MOSPD3 antibody

Specificity

Human; other species not tested.

Antibody description

MOSPD3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of MOSPD3 recombinant protein (Accession Number: XM_005250531). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112737(MOSPD3 Antibody) at dilution of 1:600