

# MTCO3 antibody

Catalog Number: 112768

## Product name

MTCO3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MTCO3 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 32-35 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: YP\_003024032). Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-COX3 (IP:Catalog No:112768, 4ug; Detection:Catalog No:112768 1:300) with mouse heart tissue lysate 4000ug.