

# MRAS antibody

Catalog Number: 112801

## Product name

MRAS antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

MRAS Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, mouse lung tissue. Positive IHC detected in human skeletal muscle. Observed molecular weight by Western-blot: 25kd

## Preparation

This antibody was obtained by immunization of MRAS recombinant protein (Accession Number: BC047690). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112801(MRAS antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:112801(MRAS antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog

# MRAS antibody



Catalog Number: 112801

---

No:112801(MRAS antibody) at dilution of 1:50 (under 40x lens)