

MRPL47 antibody



Catalog Number: 112831

Product name

MRPL47 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MRPL47 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of MRPL47 recombinant protein (Accession Number: BC032522). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112831(MRPL47 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:112831(MRPL47 Antibody) at dilution of 1:50 (under 40x lens)