

MS4A12 antibody



Catalog Number: 112864

Product name

MS4A12 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MS4A12 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of MS4A12 recombinant protein (Accession Number: NM_017716). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:112864(MS4A12 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:112864(MS4A12 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:112864(MS4A12 Antibody) at dilution of 1:50