

# MSK1 antibody

Catalog Number: 112871

## Product name

MSK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MSK1 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of MSK1 recombinant protein (Accession Number: NM\_182398). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of MCF-7 cells, using RPS6KA5 antibody Catalog No:112871 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112871(RPS6KA5 Antibody) at dilution of 1:500