

# MUC20 antibody



Catalog Number: 112890

## Product name

MUC20 antibody

## Specificity

human; other species not tested.

## Antibody description

MUC20 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human kidney tissue. Observed molecular weight by Western-blot: 70-75kd

## Preparation

This antibody was obtained by immunization of MUC20 recombinant protein (Accession Number: NM\_020790). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112890(MUC20 antibody) at dilution of 1:400