

MYCBPAP antibody



Catalog Number: 112916

Product name

MYCBPAP antibody

Specificity

Human; other species not tested.

Antibody description

MYCBPAP Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 108 kDa

Preparation

This antibody was obtained by immunization of MYCBPAP recombinant protein (Accession Number: NM_032133). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112916(MYCBPAP antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using MYCBPAP antibody Catalog No:112916 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).