

NARG1 antibody



Catalog Number: 112953

Product name

NARG1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

NARG1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 95-100 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of NARG1 recombinant protein (Accession Number: NM_057175). Purification method: Antigen Affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112953(NARG1 antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality