

MYLK antibody



Catalog Number: 112972

Product name

MYLK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYLK Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, Hela cells. Positive WB detected in mouse uterus tissue, mouse large intestine tissue. Observed molecular weight by Western-blot: 135kd

Preparation

This antibody was obtained by immunization of MYLK recombinant protein (Accession Number: BC113456). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse uterus tissue were subjected to SDS PAGE followed by western blot with Catalog No:112972(MYLK antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using MYLK antibody Catalog No:112972 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).