

MYO18A antibody



Catalog Number: 112980

Product name

MYO18A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYO18A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells, Neuro-2a cells, NIH/3T3 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 230kd, 190kd

Preparation

This antibody was obtained by immunization of MYO18A recombinant protein (Accession Number: NM_001346767). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112980(MYO18A antibody) at dilution of 1:6000



IP Result of anti-MYO18A (IP:Catalog No:112980, 3ug; Detection:Catalog No:112980 1:300) with HeLa cells lysate 3000ug.