

MYO19 antibody

Catalog Number: 112981

Product name

MYO19 antibody

Specificity

Human, Monkey; other species not tested.

Antibody description

MYO19 Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in HeLa cells, A549 cells, COS-7 cells. Observed molecular weight by Western-blot: 109 kDa

Preparation

This antibody was obtained by immunization of MYO19 recombinant protein (Accession Number: NM_025109). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112981(MYO19 antibody) at dilution of 1:1000



Immunofluorescent analysis of A549 cells using Catalog No:112981(MYO19 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)