

MYO7A antibody



Catalog Number: 112986

Product name

MYO7A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MYO7A Rabbit Polyclonal antibody. Positive WB detected in L02 cells, A431 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 160-255kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000260). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:112986(MYO7A antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using MYO7A antibody Catalog No:112986 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).