

NCF4 antibody

Catalog Number: 113035

Product name

NCF4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NCF4 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in L02 cells. Observed molecular weight by Western-blot: 41 kDa

Preparation

This antibody was obtained by immunization of NCF4 recombinant protein (Accession Number: NM_013416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113035(NCF4 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:113035(NCF4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:113035(NCF4 antibody) at dilution of 1:50 (under 40x lens)