

NCOA7 antibody



Catalog Number: 113045

Product name

NCOA7 antibody

Specificity

Human; other species not tested.

Antibody description

NCOA7 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, MCF7 cells. Observed molecular weight by Western-blot: 120-130 kDa

Preparation

This antibody was obtained by immunization of NCOA7 recombinant protein (Accession Number: XM_024446331). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113045(NCOA7 Antibody) at dilution of 1:1000



IP Result of anti-NCOA7 (IP:Catalog No:113045, 4ug; Detection:Catalog No:113045 1:1000) with HEK-293 cells lysate 3600ug.