

# NDUFB6 antibody



Catalog Number: 113075

## Product name

NDUFB6 antibody

## Immunogen

[Human NDUFB6 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NDUFB6 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 16-17 kDa

## Preparation

This antibody was obtained by immunization of NDUFB6 recombinant protein (Accession Number: NM\_002493). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113075(NDUFB6 antibody) at dilution of 1:800