

NDUFS3 antibody



Catalog Number: 113082

Product name

NDUFS3 antibody

Recommended Dilution:

WB: 1:200-1:1000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:2000


IHC: 1:20-1:200

Antibody description

NDUFS3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HEK-293 cells, HepG2 cells, human skeletal muscle tissue, mouse heart tissue, mouse skeletal muscle tissue. Positive IP detected in mouse heart tissue. Positive IHC detected in human testis tissue, human kidney tissue, human liver cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 30 kDa

IF: 1:10-1:100

Validations




Preparation

This antibody was obtained by immunization of NDUFS3 recombinant protein (Accession Number: NM_004551). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113082(NDUFS3 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human testis using Catalog No:113082(NDUFS3 antibody) at dilution of 1:50 (under 10x lens)


Applications

ELISA, WB, IHC, IF, IP


Dilutions

NDUFS3 antibody


Catalog Number: 113082



Immunohistochemical of paraffin-embedded human testis using Catalog No:113082(NDUFS3 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:113082(NDUFS3 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



IP Result of anti-NDUFS3 (IP:Catalog No:113082, 4ug; Detection:Catalog No:113082 1:500) with mouse heart tissue lysate 4000ug.