

NEK6 antibody

Catalog Number: 113105

Product name

NEK6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NEK6 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IF detected in COS-7 cells.

Preparation

This antibody was obtained by immunization of NEK6 recombinant protein (Accession Number: XM_006716936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113105(NEK6 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed COS-7 cells using Catalog No:113105(NEK6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)