

# NELF antibody

Catalog Number: 113108

## Product name

NELF antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NELF Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, human colon tissue, mouse testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 54-60 kDa

## Preparation

This antibody was obtained by immunization of NELF recombinant protein (Accession Number: BC016862). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113108(NELF antibody) at dilution of 1:100



(under 10x lens)

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113108(NELF Antibody) at dilution of 1:1000