

# NELL2 antibody

Catalog Number: 113113

## Product name

NELL2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NELL2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue. Positive IHC detected in human colon cancer tissue, human gliomas tissue. Observed molecular weight by Western-blot: 100-110 kDa

## Preparation

This antibody was obtained by immunization of NELL2 recombinant protein (Accession Number: NM\_006159). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113113(NELL2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113113(NELL2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113113(NELL2 antibody) at dilution of 1:50

# NELL2 antibody



Catalog Number: 113113

---

(under 40x lens)