

# Neogenin-Specific antibody



Catalog Number: 113114

## Product name

Neogenin-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Neogenin-Specific Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, mouse thymus tissue. Observed molecular weight by Western-blot: 160 kDa, 140 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002499). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113114(NEO1-Specific antibody) at dilution of 1:500