

NETO2 antibody



Catalog Number: 113117

Product name

NETO2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NETO2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of NETO2 recombinant protein (Accession Number: BC012381). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113117(NETO2 antibody) at dilution of 1:500