

NEURL2 antibody



Catalog Number: 113122

Product name

NEURL2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NEURL2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Positive WB detected in HL-60 cells, human heart tissue. Observed molecular weight by Western-blot: 30-32 kDa

Preparation

This antibody was obtained by immunization of NEURL2 recombinant protein (Accession Number: NM_080749). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:113122(NEURL2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:113122(NEURL2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:113122(NEURL2 Antibody) at dilution of 1:50