

NFATC3 antibody

Catalog Number: 113143

Product name

NFATC3 antibody

WB: 1:200-1:1000

IP: 1:200-1:2000

Specificity


Human, Mouse; other species not tested.

IHC: 1:20-1:200

Antibody description

NFATC3 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tissue, human testis tissue. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in mouse thymus tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 130-170 kDa

Validations




Preparation

This antibody was obtained by immunization of NFATC3 recombinant protein (Accession Number: NM_173165). Purification method: Antigen affinity purified.

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:113143(NFATC3 Antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, WB, IP

Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human ovary using Catalog No:113143(NFATC3 antibody) at dilution of 1:100 (under 40x lens)



NFATC3 antibody



Catalog Number: 113143

IP Result of anti-NFATC3 (IP:Catalog No:113143, 5ug; Detection:Catalog No:113143 1:500) with mouse skeletal muscle tissue lysate 2000ug.