

# NFIX-Specific antibody



Catalog Number: 113153

## Product name

NFIX-Specific antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NFIX-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells. Observed molecular weight by Western-blot: 23kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002501). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113153(NFIX-Specific antibody) at dilution of 1:500