

NFKB1,p105,p50-Specific antibody



Catalog Number: 113155

Product name

NFKB1,p105,p50-Specific antibody

Specificity

Human; other species not tested.

Antibody description

NFKB1,p105,p50-Specific Rabbit Polyclonal antibody. Positive WB detected in A431 cells, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 50kd,105kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003998). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113155(NFKB1,p105,p50-Specific antibody) at dilution of 1:200