

NPR2 antibody



Catalog Number: 113220

Product name

NPR2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NPR2 Rabbit Polyclonal antibody. Positive WB detected in transfected cells. Observed molecular weight by Western-blot: 111 kDa, 117 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003995). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



WB result of anti-NPR2 (Catalog No:113220, 1:1000) with HEK293 overexpressed NPR2 by Jeremy Egbert and Laurinda Jaffe. (Another anti serum 6328 as control)