

NPS antibody



Catalog Number: 113222

Product name

NPS antibody

Specificity

Human; other species not tested.

Antibody description

NPS Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human pancreas tissue.

Preparation

This antibody was obtained by immunization of NPS recombinant protein (Accession Number: NM_001030013). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113222(NPS Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain slide using Catalog No:113222(NPS Antibody) at dilution of 1:50 (under 40x lens)