

# NRIP2 antibody



Catalog Number: 113279

---

## Product name

NRIP2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NRIP2 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human kidney tissue. Observed molecular weight by Western-blot: 31kd

## Preparation

This antibody was obtained by immunization of NRIP2 recombinant protein (Accession Number: NM\_031474). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113279(NRIP2 antibody) at dilution of 1:800