

NOLC1 antibody

Catalog Number: 113299

Product name

NOLC1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Validations

Antibody description

NOLC1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human colon cancer tissue. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 130kd




Preparation

This antibody was obtained by immunization of NOLC1 recombinant protein (Accession Number: BC006769). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113299(P130 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:2000


Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113299(P130 antibody) at dilution of 1:50 (under 10x lens)



NOLC1 antibody

Catalog Number: 113299

Immunofluorescent analysis of HepG2 cells, using NOLC1 antibody Catalog No:113299 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-P130 (IP:Catalog No:113299, 4ug; Detection:Catalog No:113299 1:500) with HeLa cells lysate 2000ug.