

Oct6-Specific antibody



Catalog Number: 113319

Product name

Oct6-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Oct6-Specific Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002699). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113319(POU3F1-Specific antibody) at dilution of 1:800