

NUP98-NUP96 antibody

Catalog Number: 113386

Product name

NUP98-NUP96 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NUP98-NUP96 Rabbit Polyclonal antibody.
Positive WB detected in COLO 320 cells, HeLa cells, Jurkat cells, K-562 cells, MCF7 cells.
Positive IP detected in COLO 320 cells. Observed molecular weight by Western-blot: 98kd

Preparation

This antibody was obtained by immunization of NUP98-NUP96 recombinant protein (Accession Number: BC012906). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:113386(NUP98-NUP96 antibody) at dilution of



1:1000

IP Result of anti-NUP98-NUP96 (IP:Catalog No:113386, 4ug; Detection:Catalog No:113386 1:1000) with COLO 320 cells lysate 2000ug.