

NSUN6 antibody



Catalog Number: 113390

Product name

NSUN6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NSUN6 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells. Observed molecular weight by Western-blot: 52-55 kDa

Preparation

This antibody was obtained by immunization of NSUN6 recombinant protein (Accession Number: NM_182543). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:113390(NSUN6 Antibody) at dilution of 1:1000