

Nucleoside phosphorylase antibody



Catalog Number: 113409

Product name

Nucleoside phosphorylase antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Nucleoside phosphorylase Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, K-562 cells, Raji cells, U-937 cells. Observed molecular weight by Western-blot: 32kd

Preparation

This antibody was obtained by immunization of Nucleoside phosphorylase recombinant protein (Accession Number: NM_000270). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113409(NP antibody) at dilution of 1:500