

OSM antibody

Catalog Number: 113440

Product name

OSM antibody

Specificity

Human; other species not tested.

Antibody description

OSM Rabbit Polyclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of OSM recombinant protein (Accession Number: NM_020530). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113440(OSM Antibody) at dilution of 1:1000



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:113440(OSM Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)