

Nurr1/NR4A2 antibody



Catalog Number: 113446

Product name

Nurr1/NR4A2 antibody

Specificity

Human; other species not tested.

Antibody description

Nurr1/NR4A2 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue.

Preparation

This antibody was obtained by immunization of Nurr1/NR4A2 recombinant protein (Accession Number: XM_047444551). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113446(NR4A2 antibody) at dilution of 1:50 (under 10x lens)