

Osteopontin antibody



Catalog Number: 113507

Product name

Osteopontin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Osteopontin Rabbit Polyclonal antibody. Positive IHC detected in human stomach cancer tissue, human kidney tissue, human small intestine tissue. Positive WB detected in mouse kidney tissue, rat kidney tissue. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of Osteopontin recombinant protein (Accession Number: NM_000582). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:50-1:500


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113507(OPN, SPP1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using Catalog No:113507(OPN, SPP1 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using Catalog

Osteopontin antibody



Catalog Number: 113507

No:113507(OPN, SPP1 Antibody) at dilution of 1:200 (under 40x lens).