

# OTUD4 antibody



Catalog Number: 113514

## Product name

OTUD4 antibody

## Immunogen

[Human OTUD4 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

OTUD4 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in HeLa cells, MCF-7 cells, Raji cells. Observed molecular weight by Western-blot: 140-150 kDa

## Preparation

This antibody was obtained by immunization of OTUD4 recombinant protein (Accession Number: NM\_001102653). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113514(OTUD4 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:113514(OTUD4 Antibody) at dilution of 1:50 (under 40x lens)