

# PAX7 antibody

Catalog Number: 113608

## Product name

PAX7 antibody

## Specificity

Human; other species not tested.

## Antibody description

PAX7 Rabbit Polyclonal antibody. Positive WB detected in NCCIT cell, C2C12 cell, HeLa cells, Jurkat cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002584). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113608(PAX7 Antibody) at dilution of 1:1000



NCCIT cell were subjected to SDS PAGE followed by western blot with Catalog No:113608(PAX7 Antibody) at dilution of 1:600