

PBRM1 antibody

Catalog Number: 113613

Product name

PBRM1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PBRM1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 180-195 kDa

Preparation

This antibody was obtained by immunization of PBRM1 recombinant protein (Accession Number: BC015323). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113613(PBRM1 Antibody) at dilution of 1:1000



IP Result of anti-PBRM1 (IP:Catalog No:113613, 3ug; Detection:Catalog No:113613 1:300) with HEK-293 cells lysate 1800ug.