

PCAF antibody

Catalog Number: 113620

Product name

PCAF antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCAF Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive IP detected in HeLa cells.

Preparation

This antibody was obtained by immunization of PCAF recombinant protein (Accession Number: NM_003884). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP


Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:113620(KAT2B antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:113620(KAT2B antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-KAT2B (IP:Catalog No:113620, 4ug; Detection:Catalog No:113620 1:300) with HeLa cells lysate 2000ug.