

PDS5B antibody



Catalog Number: 113638

Product name

PDS5B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PDS5B Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells. Observed molecular weight by Western-blot: 165-170 kDa

Preparation

This antibody was obtained by immunization of PDS5B recombinant protein (Accession Number: BC039256). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:113638(PDS5B antibody) at dilution of 1:500